

APIP

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human APIP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

APIP2; CGI29; hAPIP; CGI-29; MMRP19;

Catalog #: A7102

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51074

Isotype: IgG

Swiss Prot: Q96GX9

Purity: Affinity purification

For research use only.

Background:

APIP is an APAF1 (MIM 602233)-interacting protein that acts as a negative regulator of ischemic/hypoxic injury (Cho et al., 2004)

To place an order, please [Click HERE](#).