

SEPSECS

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 55kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SEPSECS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LP; SLA; PCH2D; SLA/LP;

Catalog #: A7103

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51091

Isotype: IgG

Swiss Prot: Q9HD40

Purity: Affinity purification

For research use only.

Background:

The amino acid selenocysteine is the only amino acid that does not have its own tRNA synthetase. Instead, this amino acid is synthesized on its cognate tRNA in a three step process. The protein encoded by this gene catalyzes the third step in the process, the conversion of O-phosphoseryl-tRNA(Sec) to selenocysteinyl-tRNA(Sec).

To place an order, please [Click HERE](#).