

MTFP1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 18kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MTFP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MTP18; HSPC242;

Catalog #: A7110

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51537

Isotype: IgG

Swiss Prot: Q9UDX5

Purity: Affinity purification

For research use only.

Background:

MTP18 is a mitochondrial protein and downstream target of the phosphatidylinositol 3-kinase (see PIK3CA, MIM 171834) signaling pathway that plays a role in cell viability and mitochondrial dynamics (Tondera et al., 2004)

To place an order, please [Click HERE](#).