

METTL13

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 79kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human METTL13

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

feat; CGI-01; KIAA0859; 5630401D24Rik;

To place an order, please [Click HERE](#).

Catalog #: A7111

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51603

Isotype: IgG

Swiss Prot: Q8N6R0

Purity: Affinity purification

For research use only.