

CTNNBIP1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 9kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CTNNBIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ICAT;

Catalog #: A7122

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56998

Isotype: IgG

Swiss Prot: Q9NSA3

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene binds CTNNB1 and prevents interaction between CTNNB1 and TCF family members. The encoded protein is a negative regulator of the Wnt signaling pathway.

Two transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).