

SALL4

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 112kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SALL4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DRRS; HSAL4; ZNF797; dJ1112F19.1;

Catalog #: A7124

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57167

Isotype: IgG

Swiss Prot: Q9UJQ4

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene may be a zinc finger transcription factor. Defects in this gene are a cause of Duane-radial ray syndrome (DRRS).

To place an order, please [Click HERE](#).