

CSRNP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 64kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CSRNP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AXUD1; URAX1; TAIP-3; CSRNP-1; FAM130B;

Catalog #: A7130

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64651

Isotype: IgG

Swiss Prot: Q96S65

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that localizes to the nucleus and expression of this gene is induced in response to elevated levels of axin. The Wnt signalling pathway, which is negatively regulated by axin, is important in axis formation in early development and impaired regulation of this signalling pathway is often involved in tumors. A decreased level of expression of this gene in tumors compared to the level of expression in their corresponding normal tissues suggests that this gene product has a tumor suppressor function.

To place an order, please [Click HERE](#).