

PINK1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 63kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PINK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BRPK; PARK6

Catalog #: A7131

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 65018

Isotype: IgG

Swiss Prot: Q9BXM7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a serine/threonine protein kinase that localizes to mitochondria. It is thought to protect cells from stress-induced mitochondrial dysfunction. Mutations in this gene cause one form of autosomal recessive early-onset Parkinson disease.

To place an order, please [Click HERE](#).