

MARCKSL1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 19kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MARCKSL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

F52; MLP; MRP; MLP1; MACMARCKS;

Catalog #: A7132

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 65108

Isotype: IgG

Swiss Prot: P49006

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the myristoylated alanine-rich C-kinase substrate (MARCKS) family. Members of this family play a role in cytoskeletal regulation, protein kinase C signaling and calmodulin signaling. The encoded protein affects the formation of adherens junction. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are located on the long arm of chromosomes 6 and 10.

To place an order, please [Click HERE](#).