

ATRIP

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 86kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ATRIP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

e

Catalog #: A7139

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84126

Isotype: IgG

Swiss Prot: Q8WXE1

Purity: Affinity purification

For research use only.

Background:

This gene encodes an essential component of the DNA damage checkpoint. The encoded protein binds to single-stranded DNA coated with replication protein A. The protein also interacts with the ataxia telangiectasia and Rad3 related protein kinase, resulting in its accumulation at intranuclear foci induced by DNA damage. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).