

IL22RA2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL22RA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CRF2X; CRF2-10; CRF2-S1; IL-22BP; IL-22RA2; ZCYTOR16; IL-22R-alpha-2;

Catalog #: A7145

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 116379

Isotype: IgG

Swiss Prot: Q969J5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the class II cytokine receptor family. The encoded soluble protein specifically binds to and inhibits interleukin 22 activity by blocking the interaction of interleukin 22 with its cell surface receptor. The encoded protein may be important in the regulation of inflammatory response, and has been implicated in the regulation of tumorigenesis in the colon. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

To place an order, please [Click HERE](#).