

## ART5

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 32kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human ART5

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ARTC5;

**Catalog #:** A7146

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 116969

**Isotype:** IgG

**Swiss Prot:** Q96L15

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene belongs to the ARG-specific ADP-ribosyltransferase family. Proteins in this family regulate the function of target proteins by attaching ADP-ribose to specific amino acid residues in their target proteins. The mouse homolog lacks a glycosylphosphatidylinositol-anchor signal sequence and is predicted to be a secretory enzyme. Several transcripts encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*