

TPH2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TPH2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NTPH; ADHD7;

Catalog #: A7147

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 121278

Isotype: IgG

Swiss Prot: Q8IWU9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the pterin-dependent aromatic acid hydroxylase family. The encoded protein catalyzes the first and rate limiting step in the biosynthesis of serotonin, an important hormone and neurotransmitter. The human genome contains two related tryptophan hydroxylases, one on chromosome 11p15-p14 and one on chromosome 12q21. This gene is expressed predominantly in the brain stem. Mutations in this gene may be associated with psychiatric diseases such as bipolar affective disorder and major depression.

To place an order, please [Click HERE](#).