

PTCRA

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:29kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PTCRA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PTA; PT-ALPHA;

Catalog #:A7150

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:171558

Isotype:IgG

Swiss Prot:Q6ISU1

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a single-pass type I membrane protein that is found in immature but not mature T-cells. Along with TCRB and CD3 complex, the encoded protein forms the pre-T-cell receptor complex, which regulates early T-cell development. Four transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).