

NCR3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200 IF 1:50 - 1:100

Calculated MW: 21kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NCR3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

1C7; MALS; CD337; LY117; NKp30;

Catalog #: A7153

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 259197

Isotype: IgG

Swiss Prot: O14931

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a natural cytotoxicity receptor (NCR) that may aid NK cells in the lysis of tumor cells. The encoded protein interacts with CD3-zeta (CD247), a T-cell receptor. A single nucleotide polymorphism in the 5' untranslated region of this gene has been associated with mild malaria susceptibility. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).