

IGFL3

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 14kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IGFL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UNQ483;

Catalog #: A7154

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 388555

Isotype: IgG

Swiss Prot: Q6UXB1

Purity: Affinity purification

For research use only.

Background:

IGFL3 belongs to the insulin-like growth factor (IGF; see MIM 147440) family of signaling molecules that play critical roles in cellular energy metabolism and in growth and development, especially prenatal growth (Emtage et al., 2006)

To place an order, please [Click HERE](#).