

PDPN

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 17kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PDPN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

T1A; GP36; GP40; Gp38; OTS8; T1A2; TI1A; T1A-2; AGGRUS; HT1A-1; PA2.26;

Catalog #: A7166

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10630

Isotype: IgG

Swiss Prot: Q86YL7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a type-I integral membrane glycoprotein with diverse distribution in human tissues. The physiological function of this protein may be related to its mucin-type character. The homologous protein in other species has been described as a differentiation antigen and influenza-virus receptor. The specific function of this protein has not been determined but it has been proposed as a marker of lung injury. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).