

ADPGK

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 54kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ADPGK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADP-GK; 2610017G09Rik;

Catalog #: A7172

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83440

Isotype: IgG

Swiss Prot: Q9BRR6

Purity: Affinity purification

For research use only.

Background:

ADPGK (EC 2.7.1.147) catalyzes the ADP-dependent phosphorylation of glucose to glucose-6-phosphate and may play a role in glycolysis, possibly during ischemic conditions (Ronimus and Morgan, 2004)

To place an order, please [Click HERE](#).