

CRISP2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CRISP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

CT36; TPX1; TSP1; GAPDL5; CRISP-2;

To place an order, please [Click HERE](#).

Catalog #: A7177

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7180

Isotype: IgG

Swiss Prot: P16562

Purity: Affinity purification

For research use only.