

ELF4

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 71kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ELF4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

MEF; ELFR;

Catalog #: A7180

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2000

Isotype: IgG

Swiss Prot: Q99607

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transcriptional activator that binds and activates the promoters of the CSF2, IL3, IL8, and PRF1 genes. The encoded protein is involved in natural killer cell development and function, innate immunity, and induction of cell cycle arrest in naive CD8+ cells. Two transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).