

ELF5

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ELF5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ESE2;

Catalog #: A7181

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2001

Isotype: IgG

Swiss Prot: Q9UKW6

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of an epithelium-specific subclass of the Ets transcription factor family. In addition to its role in regulating the later stages of terminal differentiation of keratinocytes, it appears to regulate a number of epithelium-specific genes found in tissues containing glandular epithelium such as salivary gland and prostate. It has very low affinity to DNA due to its negative regulatory domain at the amino terminus. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).