

FOXN2

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FOXN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HTLF;

Catalog #: A7185

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3344

Isotype: IgG

Swiss Prot: P32314

Purity: Affinity purification

For research use only.

Background:

This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat.

To place an order, please [Click HERE](#).