

MAG

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:3000

Calculated MW: 69kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MAG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GMA; S-MAG; SIGLEC4A; SIGLEC-4A;

Catalog #: A7186

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4099

Isotype: IgG

Swiss Prot: P20916

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a type I membrane protein and member of the immunoglobulin superfamily. It is thought to be involved in the process of myelination. It is a lectin that binds to sialylated glycoconjugates and mediates certain myelin-neuron cell-cell interactions. Three alternatively spliced transcripts encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).