

GTF2H3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GTF2H3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P34; BTF2; TFB4; TFIIH;

Catalog #: A7188

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2967

Isotype: IgG

Swiss Prot: Q13889

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the TFB4 family. The encoded protein is a subunit of the core-TFIIH basal transcription factor and localizes to the nucleus. The encoded protein is involved in RNA transcription by RNA polymerase II and nucleotide excision repair and associates with the Cdk-activating kinase complex. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 14.

To place an order, please [Click HERE](#).