

IDH2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP CHIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50 IP 1:20 - 1:50

ChIP 1:20 - 1:100

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IDH2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IDH; IDP; IDHM; IDPM; ICD-M; D2HGA2; mNADP-IDH;

Catalog #: A7190

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3418

Isotype: IgG

Swiss Prot: P48735

Purity: Affinity purification

For research use only.

Background:

Isocitrate dehydrogenases catalyze the oxidative decarboxylation of isocitrate to 2-oxoglutarate. These enzymes belong to two distinct subclasses, one of which utilizes NAD(+) as the electron acceptor and the other NADP(+). Five isocitrate dehydrogenases have been reported: three NAD(+)-dependent isocitrate dehydrogenases, which localize to the mitochondrial matrix, and two NADP(+)-dependent isocitrate dehydrogenases, one of which is mitochondrial and the other predominantly cytosolic. Each NADP(+)-dependent isozyme is a homodimer. The protein encoded by this gene is the NADP(+)-dependent isocitrate dehydrogenase found in the mitochondria. It plays a role in intermediary metabolism and energy production. This protein may tightly associate or interact with the pyruvate dehydrogenase complex. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).