

KLF9

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KLF9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BTEB; BTEB1;

Catalog #: A7196

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 687

Isotype: IgG

Swiss Prot: Q13886

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transcription factor that binds to GC box elements located in the promoter. Binding of the encoded protein to a single GC box inhibits mRNA expression while binding to tandemly repeated GC box elements activates transcription.

To place an order, please [Click HERE](#).