

CDC16

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 72kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CDC16

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

APC6; CUT9; ANAPC6;

Catalog #: A7197

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8881

Isotype: IgG

Swiss Prot: Q13042

Purity: Affinity purification

For research use only.

Background:

This gene encodes a component protein of the APC complex, which is composed of eight proteins and functions as a protein ubiquitin ligase. The APC complex is a cyclin degradation system that governs exit from mitosis. Each component protein of the APC complex is highly conserved among eukaryotic organisms. This protein and two other APC complex proteins, CDC23 and CDC27, contain a tetratricopeptide repeat (TPR), a protein domain that may be involved in protein-protein interaction. Multiple alternatively spliced variants, encoding the same protein, have been identified.

To place an order, please [Click HERE](#).