

NAP1L3

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 58kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NAP1L3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MB20; NPL3;

Catalog #: A7202

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4675

Isotype: IgG

Swiss Prot: Q99457

Purity: Affinity purification

For research use only.

Background:

This gene is intronless and encodes a member of the nucleosome assembly protein (NAP) family. This gene is linked closely to a region of genes responsible for several X-linked mental retardation syndromes.

To place an order, please [Click HERE](#).