

PHIP

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 207kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PHIP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ndrp; BRWD2; WDR11; DCAF14;

Catalog #: A7207

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55023

Isotype: IgG

Swiss Prot: Q8WWQ0

Purity: Affinity purification

For research use only.

Background:

PHIP binds the pleckstrin homology (PH) domain of insulin receptor substrate-1 (IRS1; MIM 147545), modulates insulin signaling, and plays a role in pancreatic beta cell growth and survival (Farhang-Fallah et al., 2000)

To place an order, please [Click HERE](#).