

RSRC1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RSRC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BM-011; SFRS21; SRrp53;

Catalog #: A7210

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51319

Isotype: IgG

Swiss Prot: Q96IZ7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the serine and arginine rich-related protein family. The encoded protein is involved in both constitutive and alternative mRNA splicing. This gene may be associated with schizophrenia. A pseudogene of this gene is located on chromosome 9. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).