

SETBP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 175kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SETBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SEB; MRD29;

Catalog #: A7212

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26040

Isotype: IgG

Swiss Prot: Q9Y6X0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which contains a several motifs including a ski homology region and a SET-binding region in addition to three nuclear localization signals. The encoded protein has been shown to bind the SET nuclear oncogene which is involved in DNA replication. Mutations in this gene are associated with Schinzel-Giedion midface retraction syndrome. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).