

SYNCRIP

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 70kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SYNCRIP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

PP68; NSAP1; GRYRBP; HNRNPQ; HNRNPQ1; GRY-RBP; hnRNP-Q;

Catalog #: A7219

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10492

Isotype: IgG

Swiss Prot: O60506

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the cellular heterogeneous nuclear ribonucleoprotein (hnRNP) family. hnRNPs are RNA binding proteins that complex with heterogeneous nuclear RNA (hnRNA) and regulate alternative splicing, polyadenylation, and other aspects of mRNA metabolism and transport. The encoded protein plays a role in multiple aspects of mRNA maturation and is associated with several multiprotein complexes including the apoB RNA editing-complex and survival of motor neurons (SMN) complex. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the short arm of chromosome 20.

To place an order, please [Click HERE](#).