

TOB2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TOB2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TOB4; TOBL; TROB2;

Catalog #: A7223

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10766

Isotype: IgG

Swiss Prot: Q14106

Purity: Affinity purification

For research use only.

Background:

TOB2 belongs to the TOB (see TOB1; MIM 605523)/BTG1 (MIM 109580) family of antiproliferative proteins, which are involved in the regulation of cell cycle progression.

To place an order, please [Click HERE](#).