

UBE2Z

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 38kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UBE2Z

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1mg/ml

Synonym:

USE1; HOYS7;

Catalog #: A7225

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 65264

Isotype: IgG

Swiss Prot: Q9H832

Purity: Affinity purification

For research use only.

Background:

This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis.

To place an order, please [Click HERE](#).