

THRSP

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:17kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human THRSP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

S14; Lpgp; THRP; LPGP1; SPOT14;

Catalog #:A7232

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7069

Isotype:IgG

Swiss Prot:Q92748

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is similar to the gene product of S14, a rat gene whose expression is limited to liver and adipose tissue and is controlled by nutritional and hormonal factors. This gene has been shown to be expressed in liver and adipocytes, particularly in lipomatous modules. It is also found to be expressed in lipogenic breast cancers, which suggests a role in controlling tumor lipid metabolism.

To place an order, please [Click HERE](#).