

SLC35A2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 41kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SLC35A2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UGT; CDGX; UGAT; UGT1; UGT2; UGTL; CDG2M; UGALT; UDP-Gal-Tr;

Catalog #: A7233

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7355

Isotype: IgG

Swiss Prot: P78381

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the nucleotide-sugar transporter family. The encoded protein is a multi-pass membrane protein. It transports UDP-galactose from the cytosol into Golgi vesicles, where it serves as a glycosyl donor for the generation of glycans. Mutations in this gene cause congenital disorder of glycosylation type II_m (CDG2M). Multiple alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

To place an order, please [Click HERE](#).