

IGHMBP2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 109kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IGHMBP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HCSA; HMN6; CATF1; CMT2S; SMARD1; SMUBP2; ZFAND7;

Catalog #: A7240

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3508

Isotype: IgG

Swiss Prot: P40694

Purity: Affinity purification

For research use only.

Background:

This gene encodes a helicase superfamily member that binds a specific DNA sequence from the immunoglobulin mu chain switch region. Mutations in this gene lead to spinal muscle atrophy with respiratory distress type 1.

To place an order, please [Click HERE](#).