

BNIP1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 26kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BNIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NIP1; SEC20; TRG-8;

Catalog #: A7263

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 662

Isotype: IgG

Swiss Prot: Q12981

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein, which protects cells from virally-induced cell death. The encoded protein also interacts with E1B 19 kDa-like sequences of BCL2, another apoptotic protector. In addition, this protein is involved in vesicle transport into the endoplasmic reticulum. Alternative splicing of this gene results in four protein products with identical N- and C-termini.

To place an order, please [Click HERE](#).