

BLID

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 12kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BLID

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BRCC2;

Catalog #: A7275

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 414899

Isotype: IgG

Swiss Prot: Q8IZY5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a BH3-like motif containing protein involved in cell death. The encoded protein may induce apoptosis in a caspase-dependent manner. The protein is localized in both the cytoplasm and the mitochondrion.

To place an order, please [Click HERE](#).