

GCM1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GCM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GCMA; hGCMA;

Catalog #: A7282

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8521

Isotype: IgG

Swiss Prot: Q9NP62

Purity: Affinity purification

For research use only.

Background:

This gene encodes a DNA-binding protein with a gcm-motif (glial cell missing motif). The encoded protein is a homolog of the Drosophila glial cells missing gene (gcm). This protein binds to the GCM-motif (A/G)CCCGCAT, a novel sequence among known targets of DNA-binding proteins. The N-terminal DNA-binding domain confers the unique DNA-binding activity of this protein.

To place an order, please [Click HERE](#).