

GOLM1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GOLM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GP73; HEL46; GOLPH2; C9orf155; PSEC0257; bA379P1.3;

Catalog #: A7283

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51280

Isotype: IgG

Swiss Prot: Q8NBJ4

Purity: Affinity purification

For research use only.

Background:

The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. The protein encoded by this gene is a type II Golgi transmembrane protein. It processes proteins synthesized in the rough endoplasmic reticulum and assists in the transport of protein cargo through the Golgi apparatus. The expression of this gene has been observed to be upregulated in response to viral infection. Alternatively spliced transcript variants encoding the same protein have been described for this gene.

To place an order, please [Click HERE](#).