

NFIL3

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NFIL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

E4BP4; IL3BP1; NFIL3A; NF-IL3A;

Catalog #: A7294

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4783

Isotype: IgG

Swiss Prot: Q16649

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transcriptional regulator that binds as a homodimer to activating transcription factor (ATF) sites in many cellular and viral promoters. The encoded protein represses PER1 and PER2 expression and therefore plays a role in the regulation of circadian rhythm. Three transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).