

NGFRAP1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 13kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NGFRAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

Bex; BEX3; NADE; HGR74; DXS6984E;

To place an order, please [Click HERE](#).

Catalog #: A7296

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 27018

Isotype: IgG

Swiss Prot: Q00994

Purity: Affinity purification

For research use only.