

PDCD5

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 14kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PDCD5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TFAR19;

Catalog #: A7298

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9141

Isotype: IgG

Swiss Prot: O14737

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is upregulated during apoptosis where it translocates rapidly from the cytoplasm to the nucleus. The encoded protein may be an important regulator of K(lysine) acetyltransferase 5 (a protein involved in transcription, DNA damage response and cell cycle control) by inhibiting its proteasome-dependent degradation. Pseudogenes have been identified on chromosomes 5 and 12

To place an order, please [Click HERE](#).