

SETD8

Reactivity: Human Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SETD8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

SET8; KMT5A; SET07; PR-Set7;

To place an order, please [Click HERE](#).

Catalog #: A7305

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 387893

Isotype: IgG

Swiss Prot: Q9NQR1

Purity: Affinity purification

For research use only.