

SGTA

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SGTA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SGT; hSGT; alphaSGT;

Catalog #: A7306

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6449

Isotype: IgG

Swiss Prot: O43765

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which is capable of interacting with the major nonstructural protein of parvovirus H-1 and 70-kDa heat shock cognate protein; however, its function is not known. Since this transcript is expressed ubiquitously in various tissues, this protein may serve a housekeeping function.

To place an order, please [Click HERE](#).