

TFEB

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TFEB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

μg/ml

Synonym:

TCFEB; BHLHE35; ALPHATFEB;

To place an order, please [Click HERE](#).

Catalog #: A7311

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7942

Isotype: IgG

Swiss Prot: P19484

Purity: Affinity purification

For research use only.