

TWIST1

Reactivity: Human

Tested applications: WB IF

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW: 21kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TWIST1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fk

Synonym:

CRS; SCS; ACS3; CRS1; BPES2; BPES3; TWIST; bHLHa38;

Catalog #: A7314

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7291

Isotype: IgG

Swiss Prot: Q15672

Purity: Affinity purification

For research use only.

Background:

Basic helix-loop-helix (bHLH) transcription factors have been implicated in cell lineage determination and differentiation. The protein encoded by this gene is a bHLH transcription factor and shares similarity with another bHLH transcription factor, Dermo1. The strongest expression of this mRNA is in placental tissue; in adults, mesodermally derived tissues express this mRNA preferentially. Mutations in this gene have been found in patients with Saethre-Chotzen syndrome.

To place an order, please [Click HERE](#).