

TXNDC5

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 48kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TXNDC5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

ERP46; HCC-2; STRF8; PDIA15; UNQ364; ENDOPDI;

Catalog #: A7315

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81567

Isotype: IgG

Swiss Prot: Q8NBS9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein-disulfide isomerase. Its expression is induced by hypoxia and its role may be to protect hypoxic cells from apoptosis. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring upstream MUTED (muted homolog) gene.

To place an order, please [Click HERE](#).