

ETS2

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:53kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ETS2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

kq

Synonym:

ETS2IT1;

Catalog #:A7329

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2114

Isotype:IgG

Swiss Prot:P15036

Purity:Affinity purification

For research use only.

Background:

This gene encodes a transcription factor which regulates genes involved in development and apoptosis. The encoded protein is also a protooncogene and shown to be involved in regulation of telomerase. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).